## SHARPAK Infection control products





## Sharpak 220 'plus' sharps container with **VisiLOCK**

The Sharpak<sup>®</sup> 220 'PLUS' was developed following many hours of customer and market research as well as over fifteen years experience gained in the market supplying infection control products to the healthcare sector.

The SHARPAK 220 'PLUS' is approved to the British Standard BS7320:1990 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

The SHARPAK 220 'PLUS' offers the following features:

MANUFACTURING AND TESTED TO BS7320:1990

UN TESTED AND CERTIFIED

QUICK AND EASY TO ASSEMBLE

EASY TO OPERATE ROTATING OPEN/CLOSE/LOCK CLOSURE

SECURE CLOSURE LOCK

FEATURING VISILOCK INDICATOR SYSTEM

DISTINCTIVE COLOURING AND CLEAR LABELLING

ONE PIECE BASE AND HANDLE MOULDING

TOUGH AND DURABLE

SEQUENTIAL NUMBERING FOR TRACEABILITY

COST EFFECTIVE


PURPOSE DESIGNED BRACKET FEATURE FOR SIMPLE

WALL AND TROLLEY MOUNTING

The French Standard NF X 30-500 is only applicable to sharps containers with fill line capacity from 0.51 to 10 litres - therefore the Sharpak 220 PLUS is outside this standard.




## SHARPAK Infection control products









## SHARPAK<sup>®</sup> 220 'PLUS' INFORMATION

| Height                              | 383 mm        |
|-------------------------------------|---------------|
| Width                               | 297 mm        |
| Aperture                            | 85.5 x 127 mm |
| Max Capacity                        | 22 litres     |
| Capacity to fill line               | 17.6 litres   |
| SHARPAK 220 PLUS case re-order code | SU02400       |
| Quantity per case                   | 8             |



OPEN




2ND POSITION 'CLOSED'



1ST POSITION 'HALF OPEN'



3RD POSITION 'LOCKED' WITH VISILOCK INDICATOR



Supplied by:

Patron Distributors Ltd Unit 81, Baldoyle Industrial Estate Dublin 13. Ireland Tel: +353 (0)1 8391299 Fax: +353 (0)1 8391026 Email: info@patron.ie Web: www.patron.ie

