



## Sharpak Compact sharps containers with **VisiLOCK™** PATENTED

The SHARPAK® COMPACT sharps container has been specially developed to meet the growing needs of clinical staff involved in giving injections as well as for those who have to self administer medication by means of a disposable needle and syringe.

It is not always practical to have access to a safe means of disposing of used sharps and this container provides a purpose designed product which is safe, simple to use, portable and compact in size.

The SHARPAK COMPACT is approved to the British Standard BS7320:1990, French Standard NF X 30-500 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

The SHARPAK COMPACT offers the following features:

---

MANUFACTURED, TESTED AND CERTIFIED TO  
BS7320:1990 AND NF X 30-500

---

UN TESTED AND CERTIFIED

---

SUPPLIED ASSEMBLED

---

EASY TO OPERATE ROTATING OPEN/CLOSE/LOCK  
CLOSURE

---

APERTURE EQUIPPED WITH NON RETURN PETAL

---

SECURE CLOSURE LOCK

---

FEATURING VISILOCK INDICATOR SYSTEM

---

DISTINCTIVE COLOURING AND CLEAR LABELLING

---

TOUGH AND DURABLE

---

SEQUENTIAL NUMBERING FOR TRACEABILITY

---

COST EFFECTIVE

---

PURPOSE DESIGNED STAND AVAILABLE

---

ACCEPTS UP TO 5ML SYRINGES

---

SUITABLE FOR USE WITH THE SHARPAK PROCEDURE TRAY





---

ALSO AVAILABLE WITH AN OPEN APERTURE  
(NO NON-RETURN PETALS) COMPACT 'O'



There is also a demand for a sharps container with all the same safety and performance features as the SHARPAK COMPACT sharps container but where the user wishes to be more discreet. To meet this need the SHARPAK NX5 is available. The 'only' difference

between the two products is that the SHARPAK NX5 has been manufactured with a grey base and black closure and the labelling has been kept as small and as simple as possible.



**SHARPAK™**  
**COMPACT 'O'**  
with open aperture



Compact stand



OPEN



1ST POSITION 'HALF OPEN'



2ND POSITION 'CLOSED'

**SHARPAK® COMPACT INFORMATION**

Max Capacity	0.5 litres
Capacity to fill line	0.39 litres
SHARPAK COMPACT case re-order code	SU06201
SHARPAK COMPACT quantity per case	24
SHARPAK COMPACT 'O' case re-order code	SU06202
SHARPAK COMPACT 'O' quantity per case	24
SHARPAK NX5 case order no	SU06204
SHARPAK NX5 quantity per case	24
COMPACT/NX5 stand (pack quantity 1)	SM07500

Supplied by:

Patron Distributors Ltd  
Unit 81, Baldoyle Industrial Estate  
Dublin 13, Ireland  
Tel: +353 (0)1 8391299  
Fax: +353 (0)1 8391026  
Email: info@patron.ie  
Web: www.patron.ie

