

# HOVERJACK®

Air Patient Lift

## HoverJack® Air Patient Lift

Recognized as the industry standard for safe patient lifting, the HoverJack® air patient lift allows caregivers to safely lift patients who have fallen without gathering a lift team. The HoverJack® inflates to lift patients from the floor to bed or stretcher height in a supine position, maximizing patient comfort and minimizing the risk of injury to the patient and their caregivers. Once the HoverJack® is inflated, the HoverMatt® air transfer mattress can be used with the same air supply unit for safe and easy lateral patient transfer to a bed or stretcher.



4th Chamber Inflated



3rd Chamber Inflated



2nd Chamber Inflated



1st Chamber Inflated



### HoverJack® Air Patient Lift Benefits:


- Minimize staff members needed to safely lift patients who have fallen
- Control height through individual chamber inflation
- Self-sealing inflation valves for ease of use
- Quick release deflation valves to lower patient if needed
- Teflon® bottom allows for easy movement across all surfaces
- Available in two widths with a weight capacity of 1200 lbs. to accommodate most patients



# HOVERJACK®

Air Patient Lift

HoverJack® Air Patient Lift



## Product Specifications

### Construction

Topside Material: Nylon Oxford

Underside Material:

Teflon® impregnated polyester

Latex Free

### Dimensions

Two Sizes Available:

- 32" W x 78" L – HJ3202
- 39" W x 78" L – HJ3902

### Inflation

- Four inflation chambers with self-sealing valves
- Four quick release deflation valves (red capped)

### Compatible with HoverTech Air Supplies

- AIR200G (800W), AIR400G (1100W) and HTAIR 1200 (Variable Speed)

For more information, contact your local HoverTech representative or reach us at:

800.471.2776

[www.HoverMatt.com](http://www.HoverMatt.com)

[info@hovermatt.com](mailto:info@hovermatt.com)

Ask about HELPS<sup>SM</sup> Healthcare Ergonomic Lifting Program, HoverTech's safe patient handling program that offers guaranteed results.



513 S. Clewell Street  
Bethlehem, PA 18015

**Phone** 800.471.2776

**Fax** 610.694.9601

[www.HoverMatt.com](http://www.HoverMatt.com)