

FERNO DECON

Chair & Stretchers

Offers maximum comfort during the decontamination process



FOR MORE DETAILED
INFORMATION

Specifications

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Decon Equipment** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

	Length (folded)	Width (folded)	Depth (folded)	Height (folded)	Weight	Load Limit
Model 11-D	1800mm(940mm)	600mm(600mm)	430 mm	100 mm	10.5 kg	159kg/25st/350lbs
Model 108-T	2340mm(1180mm)	585mm(200mm)	N/A	240mm(180mm)	8 kg	544kg/85st/1119lbs
Decon Chair	N/A	435mm(435mm)	430mm(235mm)	935mm(775mm)	9.1kg	200kg/31st/440lbs

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**
United Kingdom

FERNO **DECON** Chair & Stretchers

Model 11 Decon Stretcher

Offers maximum comfort during the decontamination process

This heavy duty, folding stretcher is as easy to handle as it is durable. Made of lightweight aluminium with a vinyl-coated nylon cover, it features an exclusive slimline centre hinge system. This allows it to be easily folded away for compact storage.




Model 108-Tactical Litter

Patient adopts a comfortable seated position during decontamination

The Ferno MIR Tactical Litter 108-T has strong but lightweight aluminium poles, hinges, legs and spreader bars. Handles have a non-slip surface design and are an integral part of the 6061 tubular aluminium litter poles for greater strength.

- Solid sanoprene rubber end plugs keep any blood, body fluids, water or hazardous chemicals from entering the litter poles.
- Extra strong 2024-T4 aluminium for greater strength and durability.
- Built to US Military Spec MIL-L-49511B.



Decon Chair

Patient adopts a comfortable seated position during decontamination

Specially designed for the comfortable evacuation of patients who may be contaminated as well as injured, this proven chair is ideal for transporting patients away from various hazardous environments.

- Covered in a mesh material specifically designed for decontamination
- One-piece back and seat panel allows the patient to stay seated during the decontamination process
- Durable and robust – made from stainless steel



Part Nos

015317200	Model 11 Decon Stretcher (2 wheels/2 posts)
010797000	Model 108-T Tactical Litter
071481301/D	Compact 1-Decon Chair

FERNO - Decon Chair & Stretchers FPC001 10.14

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.

 **FERNO**
United Kingdom